



Sleep Therapy



Yes, you can.

Twilight™ II



Six (6) forehead positions provide maximum comfort.



Resized breathable headgear with soft hook and loop fasteners allows for maximum adjustability and comfort.



Quick-release clips make it easy to remove the mask without changing headgear settings.

Twilight™ II' features for your sleep lab:

- Multi-patient use for cost effective care
- Easy disassembly and cleaning
- Built in port for convenient measurements

Twilight II features for your patients:

- Easy replacement cushion, available in 3 sizes
- Breathable headgear for increased comfort
- Six-position forehead support for a customized fit
- Extremely lightweight
- 360 degree ultra-quiet main swivel



Yes, you can.®

Invacare® Twilight™ II Mask “GREAT FIT” Guarantee

Invacare is proud to offer an unconditional, 30-Day No Risk “Great Fit Guarantee” to all providers of the Invacare® Twilight™ II Mask. If for any reason a patient is dissatisfied within the first 30 days of owning their Twilight II mask, Invacare will gladly send the provider a replacement mask at no charge!

Independent testing shows that patients prefer the Twilight II Mask over their usual mask.¹ We are confident that the Twilight II is a superior mask that will fit the majority of CPAP patients.

The program is as easy as 1, 2, 3:

- STEP 1:** In the event that the patient chooses to return a Twilight II Mask to the provider, they may do so within 30 days of purchase.
- STEP 2:** The provider must call the Invacare Customer Service Department at (800) 333-6900 to generate an RA and return the mask.
- STEP 3:** Invacare will gladly send a new Twilight II Mask to the provider upon receipt of the original mask. In this way, the provider keeps a Twilight II Mask in their inventory to fit another patient².



MODEL SELECTION GUIDE

ISP2100STD	Invacare® Twilight™ II mask shell, standard size nasal cushion, 4 clips and headgear
ISP2100SM	Invacare Twilight II mask shell, small size nasal cushion, 4 clips and headgear
ISP2100L	Invacare Twilight II mask shell, large size nasal cushion, 4 clips and headgear
ISP2110STD	Replacement standard size nasal cushion for Twilight II mask
ISP2110SM	Replacement small size nasal cushion for Twilight II mask
ISP2110L	Replacement large size nasal cushion for Twilight II mask
1122733	Replacement headgear and 2 clips for Twilight II mask
1122734	Chin Support for Twilight II mask

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Invacare, the Blue Medallion Design and Yes, you can. are registered trademarks and Twilight is a trademark of Invacare Corporation.

Specifications are subject to change without notification.

© 2006

Invacare Corporation
Form No. 06-030
Rev. 08/06

1. “Patient Perception of a New Nasal Mask – Results of a Clinical Trial”, Published in SLEEP, 2003; 26 (Abstract Supplement): A1213

2. Limit of one replacement mask per patient.