



- Full Efficacy Data Reporting
- Innovative tracking technology ensures accurate and comfortable therapy
- 4-25 cmH₃O therapy pressure
- Real time audible alert for tubing/mask leak
- REslex® exhalation relief with four comfort settings
- Remote compliance data access with iCode[®]
- SD Card standard
- Integration with cloud based iCode Connect[™]
- · Embedded memory stores last night's full raw data, 365 nights of user data and more than 30 years of therapy records
- Last night's compliance-therapy snapshot to assist patient with compliance
- Automatic leak and altitude compensation
- Inspiration trigger for automatic start-up
- Quiet operation

Ergonomic Design

- Integrated heated humidifier
- DC 24V powered and infrared controlled humidifier for patient safety and comfort
- Delay-off feature protects the device from humidity hazard
- Backlit LCD display for operation in dark
- Light user buttons with lockup capability



B7000 Bundled RESmart® BPAP + Humidifier B5500 RESmart® BPAP Integrated Humidifier M1020 RESmart® BPAP Air Filter (2 Pack) M1019B RESmart® Air Filter Cover - BPAP M4510 Humidifier Water Chamber

Dimensions: $8.66 \times 7.6 \times 4.4''$ (220 x 194 x 112mm), 12.3" (313 x 194 x 112mm) (with Integrated Heated Humidifier)

Weight: __ lbs. (2.2 kg), <4.8 lbs. (3 kg) (with Integrated Heated Humidifier)

Environmental

Operating Transport and Storage 5 to 35°C -20 to 55°C Temperature ≤ 80% Relative ≤ 93% Humidity Non-condensing Non-condensing

Atmospheric

860 to 1060 hPa 500 to 1060 hPa Pressure

Standards Compliance:

IEC 60601-1 General Requirements for Safety of Medical **Electrical Equipment**

IEC 60601-1-2 Electromagnetic Compatibility ISO 17510 Sleep Apnea Breathing Therapy

ISO 8185 General Requirements for Humidification Systems

Electrical:

AC Power Consumption: 100 - 240V AC, 50/60Hz, Max 95VA Type of Protection Against Electric Shock: Class II Equipment Degree of Protection Against Electric Shock: Type BF Applied

Degree of Protection Against Ingress of Water: Drip Proof,

Mode of Operation: Continuous

RESmart® Pressure **Function Auto Pressure** BPAP 25 4-25cmH₃O S, CPAP No BPAP 25A 4-25cmH₂O S, CPAP, AutoS available

Pressure

IPAP 4 to 25 cmH₃O* ±1 cmH₂O** **EPAP** ±1 cmH₂O** 4 to 25 cmH₂O* ±1 cmH₂O** CPAP 4 to 20 cmH₂O Ramp Duration 0 to 60 minutes ±10% of the setting Rise Time 0 to 3** +25%**

* Limited to 25 cmH₂O when using the Auto feature (only available under

** Dynamic pressure accuracy is ± 1 cmH2O measured at the patient end of the circuit.

Static pressure accuracy is ± 0.5 cmH2O measured at the patient end of the *** The range of values corresponds to tenths of seconds (e.g., a setting of 4

indicates a Rise Time of 0.4 seconds).

**** Measured at the patient end of circuit.

Disposal: Dispose of the device in accordance with local regulations.

Manufacturer: BMC Medical Co., Ltd.

NOTE: specifications subject to change without prior notice.

3B Medical, Inc. www.3Bproducts.com